




## Safety Data Sheet

### Section 1: Identification

Product Name: Latex Medical Exam Gloves  
Product Number: LG100 Series  
Date Prepared: 1/1/2019  
Distributor: 

### Section 2: Hazards Identification

### Section 3: Composition/Information on Ingredients

Ingredient	%
Natural Rubber Latex	≥95.5
Sulphur	<2.3
ZnO	<2.1
ZDBC	<1.5
ZDEC	<1.5
Ammonia Solution	<0.9
TiO <sub>2</sub>	<0.8
Wingstay L Powder	<0.7

### Section 4: First-Aid Measures

Inhalation: No need for first aid is anticipated  
Eye Contact: No need for first aid is anticipated  
Skin Contact: No need for first aid is anticipated  
Ingestion: No need for first aid is anticipated

### Section 5: Fire-Fighting Measures

**Suitable (and unsuitable) Extinguishing Media:** Ordinary combustible material (e.g., water, foam, carbon dioxide).

## Section 6: Accidental Release Measures

## Section 7: Handling and Storage

**Precautions for Safe Handling:** This product is considered to be an article which does not release or otherwise result in exposure to a hazardous chemical under normal use conditions.

**Conditions for Safe Storage:** Store in cool and dry place. Avoid excessive heat.

## Section 8: Exposure Controls / Personal Protection

Not Applicable

## Section 9: Physical and Chemical Properties

Form:	Elastic elastomer
Color:	Inherent color
Odor:	Inherent odor
Density:	0.92 g/cm <sup>3</sup>
Solubility:	Insoluble
Soluble In:	Some petroleum base solvents
pH:	Neutral

## Section 10: Stability and Reactivity

Stability: Stable

Materials and Conditions to Avoid: Petroleum base solvents (Sulphuric Acid 98%, Diethyl phthalate, Methacrylate)

Conditions to Avoid: Sun light, UV light, and heat.

Hazardous Decompositions Products: SO<sub>2</sub>, CO<sub>2</sub>

## Section 11: Toxicological Information

Skin Irritation: No irritation

Dermal Contact Sensitization: No sensitization

## Section 12: Ecological Information

## Section 13: Disposal Considerations

Dispose in accordance with all state, local, and federal regulations.

## Section 14: Transport Information

No special requirement

## Section 15: Regulatory Information

Regulation:	FDA CFR Title 21
Classification:	Class I non-sterile medical device
Risk Phrases:	<i>Not Applicable</i>

## Section 16: Other Information

This Product Safety Data Sheet is offered solely for your information. 3D Dental provides no warranties, either express or implied, concerning the safe use of this product in your process or in combination with other substances and assumes no responsibility for the accuracy or completeness of the data contained herein. User has the sole responsibility to determine the suitability of the product for any use and the manner of use contemplated