

VISIONARY

PSP PHOSPHOR STORAGE PLATES

- 3D Dental Phosphor plates are a digital imaging system that reads the X-ray on a phosphor imaging plate and send the image to your computer.
- These wireless plates are as thin as a business card.
- With the proper care and cleaning, the plates can be used 500-700 times before needing replacement.
- The plates feature rich images of diagnostic quality and are available in seconds.
- Choose an economical and easy to use option.



Phosphor Imaging Plates Compatible with Air Techniques

Size #0, Box of 2..... #AT-0
Size #1, Box of 2..... #AT-1
Size #2, Box of 4..... #AT-2
Size #3, Box of 2..... #AT-3
Size #4, Box of 2..... #AT-4
Size #5 x 12, Box of 1... #AT-512

Phosphor Imaging Plates Compatible with Gendex

Size #0, Box of 2..... #GT-0
Size #1, Box of 2..... #GT-1
Size #2, Box of 4..... #GT-2
Size #3, Box of 2..... #GT-3
Size #4, Box of 2..... #GT-4

Phosphor Imaging Plates Compatible with Soredex

Size #0, Box of 2..... #ST-0
Size #1, Box of 2..... #ST-1
Size #2, Box of 4..... #ST-2