


BIOVIVA

SILK SUTURES



Bioviva Silk Black Braided Sutures exemplify the quality and reliability you seek in a surgical product. They feature super fine needles as well as high knot comfort and security. Only AISI 300 Stainless steel needles are used for enhanced needle integrity, sharpness, insertion ease and patient comfort. Bioviva Silk has a two year tensile strength retention rating for maximum peace of mind.

3D-SBB1679

3-0 Reverse Cutting 19mm 3/8" 45cm
Compare to Look 784B and Ethicon 1679H

3D-SBB632

3-0 Reverse Cutting 22mm 1/2" 45cm
Compare to Look 785B and Ethicon 632G

3D-SBB684

3-0 Reverse Cutting 24mm 3/8" 45cm
Compare to Look 786B and Ethicon 684H

3D-SBB683

4-0 Reverse Cutting 19mm 3/8" 45cm
Compare to Look 781B and Ethicon 683G

dream



bioviva

CANDI-CAINE

visionary



pristine



Hydrolux

KANGAROO



THE ESSENTIALS



SABUR