


bioviva HEMOSTATIC gauze



- Innovative soluble, sterile material made from regenerated cellulose
- For topical use; Ideal for extraction sites and surgical procedures
- Upon contact with blood or fluid rapidly transforms into viscous gel, sealing wound voids and promoting clotting
- Efficient and cost effective
- User-friendly, stays in place after packing
- 100% natural. Derived from plants



Puncture foil backing and remove one piece of pre-cut fabric



Place BioViva onto the wound or extraction site with forceps



Within 15 seconds, BioViva becomes viscous

#BHG
20/Box