


SADDLE STOOL. DOCTOR STOOL. DENTAL STOOL. MEDICAL STOOL.



Back rest:

- PU Foaming(PVC)
- Two-way adjustable backrest. Allows vertical alignment for proper lumbar support and dynamic horizontal adjustment
- Back rest adjustable controls/lock by lever switch underneath the seating

Seating (Upholstery):

- PU Foaming(PVC)
- High density foam ensures shape of the seat stays unchanged for long term usage
- High Antimicrobial/bacterial and mildew resistant
- Superior abrasion resistance
- Height adjustable
- Tilt function available (3-way lever switch)

Gas lift:

- Base: Aluminum
- Gas Lift: 140mm travel

Order #CSS

CADENCE

by

3D  Dental