

UP TEMP

TEMPORARY CEMENT

- Available in both Eugenol and Non-Eugenol forms.
- The excellent dispensing system eliminates the need for mixing and any possible proportion errors.
- We love the working time and amazing handling characteristics.
- The film thickness is below 15 microns and testing has shown over 7,000 PSI in compressive strength.
- This is a great solution for temporaries where the crown is being bonded on the tooth.
- It mixes easily with no mess as well as sets and cleans beautifully.

Box contains syringe and 6 mixing tips



dream



bioviva

CANDI CAINE

visionary



pristine



Hydrolux

KANGAROO



THE ESSENTIALS



SABUR