


bioviva

HEMOSTATIC

gauze



- Innovative soluble, sterile material made from regenerated cellulose
- Upon contact with blood or fluid, rapidly transforms into viscous gel, sealing wound voids and promoting clotting
- Efficient and cost effective
- User-friendly, stays in place after packing
- 100% natural. Derived from plants



Puncture foil backing and remove one piece of pre-cut backing



Place BioViva onto the wound or extraction site with forceps



Within 15 seconds, BioViva becomes viscous

#BHG
20/Box

dream



bioviva

CANDI-CAINE

visionary



pristine



Hydrolux

KANGAROO



THE ESSENTIALS



SABUR